

16S A Site rRNA

5'-GGCGUCA CACC U
 CCCUGA GUGG C
 A G

Control RNA

5'-GGCGUCACACC U
 CCGCAG UGUGG C
 A G

Figure 1. Sequence and structure of 27mer RNA target

Patent No.: IBIS-0369
MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 2 of 33




Figure 2. MS/MS of control RNA/DNA (upper); control+paromomycin (lower)

Docket No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and Steven Hofstadler
Atty: Paul K. Leggaard - Telephone: 215 568 3100
Sheet 3 of 33




Figure 3. RNA/DNA chimera+paramomycin (upper); chimera+library (lower)

Docket No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 4 of 33




Figure 4. MS-MS analysis of member bound to RNA/DNA chimera

Docket No.: IBIS-0369
 Title: MASS SPECTROMETRIC METHODS FOR
 BIOMOLECULAR SCREENING
 Inventors: Stanley T. Crooke, Richard Griffey and
 Steven Hofstadel
 Atty: Paul K. Legaard - Telephone: 215 568 3100
 Sheet 5 of 33




Figure 5. ESI-MS of RNA/DNA chimera bound to paromomycin and library

cket No.: IBIS-0369
e: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 6 of 33




Figure 6. MS/MS of RNA/DNA chimera + compound with mass 665.1 not bound at the A-site

Doc. No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 7 of 33




Figure 7. MS-MS analysis of member bound to RNA/DNA chimera at the A-Site




Figure 8. MS Fragmentation of DNA:DNA duplexes





Figure 9. MS Fragmentation of DNA:RNA duplexes

**MASS Analysis of Binding Location
non-A site binder**


FIGURE 10



Patent No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 11 of 33


MASS Analysis of Binding Location non-A site binder

FIGURE 11



Docket No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 12 of 33

**MASS analysis of 16S A site RNA plus
216 member library
(performed on quadrupole ion trap)**



High Precision ESI-FTICR Mass Measurement of 16S A site RNA/Paromomycin Complex

use of unbound RNA as internal mass standard
provides low ppm mass measurement errors





Figure 13

Figure 14
MASS of 60-Member Ibis Library Against 16S A-site RNA



ocket No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 15 of 33

Figure 15
MASS of 60-member Library against 16S A-site Model




Figure 16
FT-ICR MS of Starting Library




Figure 17
**Compound Identification from a 60-member
Combinatorial library with MASS**

Complex M_{meas}	$9320.300 \pm .009$ Da
RNA M_{meas}	<u>$8922.189 \pm .009$</u>
ΔM	<u>$398.111 \pm .009$ Da</u>

$C_{15}H_{16}N_4O_2F_6$ 398.117 Da
 $-C_{19}H_{22}N_6O_4$ - - - - - 398.163 -

Elemental Composition Constraints

"Scaffold"
 $C_{12}H_{23}N_3O_7$


Measured Mass: 615.2969

Mass Tolerance: 1.0 ppm

Charge: 0

Element Min. atoms Max. atoms

^{12}C	12	30
1H	23	60
^{16}O	7	20
^{14}N	3	20




Possible Elemental Compositions:

Calc. Mass (amu)	Error (ppm)	Molecular Formula
615.296291	0.98	$^{16}O_4 \ ^{14}N_1 \ ^{12}C_7 \ ^1H_{33}$
615.296298	0.98	$^{16}O_9 \ ^{14}N_1 \ ^{12}C_2 \ ^1H_{39}$
615.296305	0.97	$^{16}O_{14} \ ^{14}N_5 \ ^{12}C_3 \ ^1H_{45}$
615.296808	0.15	$^{16}O_{15} \ ^{14}N_7 \ ^{12}C_8 \ ^1H_{41}$
615.296815	0.14	$^{16}O_{20} \ ^{14}N_1 \ ^{12}C_9 \ ^1H_{47}$


Further constrain by
 elemental
 composition of
 "letters"
 unintended
 products...

Figure 19
MASS K_d determination for 16S-Paromomycin




Multitarget Affinity/Specificity Screening

Figure 20




MASS Analysis of 27 Member Library with 16S A-site RNA

Figure 21.



MASS Protection Assay

Figure 22



MASS Protection Assay

Figure 23





Figure 24
Eukaryotic and Prokaryotic A-Sites
Aminoglycoside antibiotics bind to
A-site of decoding region in 16S RNA



$\Delta MW = 15.011 \text{ Da}$

Figure 25
Neutral Mass Tag Does Not Affect Ligand Binding





Figure 26
**Simultaneous Screening of 16S A-site and 18S A-site
Model RNAs Against Aminoglycoside Mixture**



Docket No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 27 of 33

Figure 27

5'

G-C
G-C
C-G
G-C
U-U
C-G

¹⁴⁰³A-A


C-G ^{A₁₄₀₃}

A-U
C-G
C-G

U U G
 C C

R A = adenosine
C A = deoxyadenosine

Figure 28



Docket No.: IBIS-0369
Title: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstadler
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 29 of 33

Figure 29




Figure 30





Figure 31



16628-1-3 16628-1-4 16628-1-5 16628-1-6

Figure 32

cket No.: IBIS-0369
: MASS SPECTROMETRIC METHODS FOR
BIOMOLECULAR SCREENING
Inventors: Stanley T. Crooke, Richard Griffey and
Steven Hofstader
Atty: Paul K. Legaard - Telephone: 215 568 3100
Sheet 32 of 33

nr	name	apex	start	stop	height	area
1	16628-1-4	1783.710	1783.635	1783.834	14.55	1.63
2	16628-1-3	1783.909	1783.834	1783.972	60.04	5.15
3	16628-1-2	1784.109	1784.021	1784.184	115.60	11.14
4	16628-1-1	1784.308	1784.233	1784.383	167.34	15.89
5	16628-1-0	1784.508	1784.433	1784.620	133.94	14.74
6	16628-1-1	1784.707	1784.620	1784.795	136.60	13.38
7	16628-1-2	1784.907	1784.795	1784.982	82.63	8.56
8	16628-1-3	1785.107	1785.032	1785.219	57.81	5.21
9	16628-1-4	1785.306	1785.232	1785.369	32.31	2.65
10	16628-1-5	1785.506	1785.456	1785.569	17.67	1.12
11	16628-10019-4	1906.974	1906.874	1907.031	12.63	1.00
12	16628-10019-3	1907.173	1907.045	1907.273	22.54	2.11
13	16628-10019-2	1907.373	1907.287	1907.444	33.86	2.91
14	16628-10019-1	1907.572	1907.458	1907.701	34.87	3.30
15	16628-10019-0	1907.772	1907.701	1907.843	20.93	1.55
16	16628-10019-1	1907.972	1907.900	1908.043	21.03	1.55
17	16628-10019-2	1908.157	1908.086	1908.271	10.97	0.90
18	16628-4	2229.874	2229.679	2230.029	27.51	4.87
19	16628-3	2230.146	2230.029	2230.263	111.72	16.23
20	16628-2	2230.380	2230.263	2230.516	225.18	32.39
21	16628-1	2230.633	2230.516	2230.770	280.66	40.90
22	16628-0	2230.887	2230.770	2231.023	287.24	41.95
23	16628-1	2231.140	2231.023	2231.257	242.23	34.17

Graph Table

Figure 33

